

## Wheelchair to IV Pole Quick Release Clamping System

**Product Code: CH-1003QR** (Patent Pending)

- One “Clamping Arm” for multiple wheelchairs
- Easy to install and remove
- Decreases wrist, arm and back strain
- “Anchor” attaches to either handle
- 9 Adjustable locking positions
- 240 degree adjustable radius
- Safe and secure



- Decreases injuries to nurses, transporting personnel and patients
- Accepts standard IV Pole diameters between ¾” to 1 ¼”



The **CH-1003QR: Wheelchair to IV Pole Quick Release Clamping System** attaches to your existing wheelchairs and IV Poles.



Attach the “Anchors” to the handle bars of all of your wheelchairs.

Once the “Anchors” have been installed; the “Clamping Arm” can be attached and removed as needed.



Simply adjust the “Clamping Arm” around the “Anchor”.



Slide the “Pin” through the “Anchor” and secure the clip.



The **CH-1003QR** attached to a standard IV Pole.

