

INTRODUCING OUR NEWEST SAFE PATIENT HANDLING CLAMP




Developed to Save: PATIENTS ♦ STAFF ♦ TIME ♦ MONEY

aXtraHand LLC
51528 Industrial Drive
New Baltimore, MI 48047
P – 586-777-0018
F – 866-777-0535
www.axtrahand.com

The Universal Clamp to IV Pole

Product Code: **CH-6000** (Patent Pending)

- ✓ 9 Adjustable locking positions
- ✓ 240 degree adjustable radius
- ✓ Accepts standard IV Pole diameters between 3/4" to 1-1/4"
- ✓ Can be easily turned 90° by removing one screw for mounting horizontally or vertically using a 5/32 allen wrench
- ✓ Attaches to any type of round or square bar on most patient transport equipment (from 1" to 2")



The **CH-6000 Universal Clamp to IV Pole HDB Clamp** wraps around the entire range of square or round bars, from 2" OD down to 1" OD, all without marring the equipment. The unique 2 position latch on the link saves time when switching from a 2" bar to a small bar, just connect to the second latch for the smaller sizes. The quick release handle takes only 2 turns to be able to drop the locking mechanism, and release the clamp from the bar.

The **IV Pole Clamp** locks in and releases IV Poles easily and quickly by a simple twist of the knob.

Make sure the handle is snug and secured tightly in the notched in area



Attached to a Square / Horizontal bar

Attached to a Round / Vertical bar